

ATMOS® Hygrotherm plus

# When precision counts!



## ATMOS® Hygrotherm plus

ATMOS® Hygrotherm plus ear irrigation		REF
1	<p><b>Irrigation unit for the ENT field</b> with electronically controlled flow heater, irrigation mode (37°C, up to 450 ml/min), energy-saving mode to reduce energy consumption, highest temperature accuracy and safety for the patient. Water supply hose G3/4", 3 m, waste water hose G1/4", 3 m, 220-240 V~, 50/60 Hz. <i>Included in delivery:</i> 3 x nozzle attachment -II straight, 80 mm; 1 x Plexiglas splash protection; 1 x seal for splash protection; 1 package of hose tips for nozzle (30 pcs.); 1 x water filter housing R.V.-A1/4"; 1 x water filter CellPor-Pe; 1 x mains cable.</p>	502.0901.0
	ATMOS® Hygrotherm plus (special voltage 100-127 V~)	502.0901.1
	<p><b>ATMOS® Aqua clean for rinsing and stimulation devices</b> pressure-free water separation in accordance with EN 1717, water decontamination and discalcification <i>Included in delivery:</i> Device, connection material, 1 litre of ATMOS® Green &amp; Clean WK disinfectant</p>	502.1200.0
Accessories		
2	<p><b>Water prefilter for ENT units, can be flushed out</b> <i>Included in delivery:</i> 1 x reduction flange G3/4"; 1 x nut - lip R1" T-chrome; 1 x solder fitting 90°; 1 x bracket for water filter; 2 x seal RD 17x24x2; 2 x seal RD 21x30x2.</p>	502.0990.0
3	Rinsing connection for thermic stimulation, teflon application help for ear irrigation	501.0355.0
4	Jet connection, for water irrigation handle, L = 80 mm	508.0427.0
	Jet connection, for water irrigation handle, L = 110 mm	508.0429.0
5	Splash protection, plexiglass	501.0331.0
	Water supply hose G3/4", 1.5 m extension	501.0315.6
	Waste water hose R 1/4", special length, surcharge/0.5 m	502.0881.1
Consumables		
	Hose tips, 30 pcs. (single use), minimum qty.: 2 packs (price per pack)	502.0844.0
	O-ring for splash protection	501.0331.1
	Filter cartridge for water filter	502.0891.0
6	<p>Hygiene filter for water flushing End positioned membrane filter for connection to the irrigation handle. Medicinal product acc. to 93/42/EWG. Service life 4 weeks (Luer Lock).</p>	000.0918.0



1 ATMOS® Hygrotherm plus




2 Water prefilter for ENT units



3 Rinsing connection for thermic stimulation



4 Jet connection for water irrigation handle



5 Splash protection



6 Hygiene filter for water flushing



**MedizinTechnik**

ATMOS MedizinTechnik GmbH & Co. KG  
Ludwig-Kegel-Str. 16  
79853 Lenzkirch / Germany  
Tel: +49 7653 689-370  
atmos@atmosmed.de

[www.atmosmed.com](http://www.atmosmed.com)