

Quality in detail

ATMOS[®] HD Endoscopes

Provide the perfect insight



Coordinated components

The quality of endoscopic systems depends largely on the perfect coordination of its components. With their high-quality optics and durable design, ATMOS HD endoscopes are the perfect addition to your high-resolution camera system:

- HD optics - for razor-sharp images
- Homogeneous illumination - from the centre to the edge of the image
- Autoclavable - for safe and efficient processing
- High-quality design - stainless steel and sapphire glass increase longevity
- Compatibility - can be used with all common endocopy cameras and light guides

Rigid endoscopes

	REF	α	\emptyset	\longleftrightarrow
Laryngoscope HD	950.0254.0	70°	8 mm	L= 190 mm
	950.0255.0	70°	10 mm	L= 190 mm
	950.0256.0	90°	10 mm	L= 190 mm
Otoscope HD	950.0257.0	0°	4 mm	L= 50 mm
	950.0258.0	30°	4 mm	L= 50 mm
	950.0259.0	0°	2,7 mm	L= 50 mm
Nasopharyngeal endoscope HD 4 mm (Sinuscope)	950.0260.0	0°	4 mm	L= 175 mm
	950.0261.0	30°	4 mm	L= 175 mm
	950.0262.0	45°	4 mm	L= 175 mm
	950.0263.0	70°	4 mm	L= 175 mm
Nasopharyngeal endoscope HD 2,7 mm (Sinuscope)	950.0264.0	0°	2,7 mm	L= 110 mm
	950.0265.0	30°	2,7 mm	L= 110 mm

Flexible endoscopes

Nasopharyngoscope	950.0243.0	0°	3,2 mm	L= 300 mm
• 80° angle of view				
• Angle 125°/125°				
• 16.000 Pixel				

Accessories

light fibre optic cable	straight	950.0152.0	0°	4,8 mm	L = 1,8 m
Storz-connection	straight	508.0663.0	0°	3,5 mm	L = 1,7 m
	90°	508.0664.0	90°	3,5 mm	L = 1,8 m